TITANIUM BINDING SYSTEM

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Tension Band Fixation Technology

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江苏双羊医疗器械有限公司
Jiangsu Shuangyang
Medical Instrument Co., Ltd.

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Titanium Binding System Introduction

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Solid and liquid all have surface tension to against fracture. So, titanium cable will have better static strength and fatigue strength along with the increase of strands.

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Strong, durable and soft. All these characteristics give Titanium Binding System great performance to be ideally taut and fixed, to make the fragments remain stable. Completely avoid the loop or kink as hard steel wire.

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
Based on the anatomical and functional purpose, tension band fixation technology of Titanium Binding System has been clinically applied to: patella fractures, olecranon fractures, proximal and distal ulna fractures, periprosthetic fractures, humerus and ankle fractures, medial malleolus fracture, pilon fracture, acromioclavicular dislocation...etc. All these fractures are characterized by obvious fracture displacement and dysfunction. Treatments of these fractures request balancing out muscle strength, but the fragments are too small to be fixed by big internal implants. So, titanium cable can play an irreplaceable role.

Titanium binding system can play an important role in many other case, such as PFF, comminuted fracture of femoral shaft, nonunion owing to failed internal fixation, reconstruction of bone defect and wide-bound splitting fracture. If need other measures to fix, titanium binding system can coordinate regular internal fixation to get better stability.

柔韧 Flexible 稳定 Stable

带孔可折断
bone needle with hole (can be broken)

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Titanium cable is made of multi-strand wire twist together, with good flexibility


产品编号 Item No.	规格 Specification	
4232	Φ 1.3	600mm
4234	Φ 1.8	600mm

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Different size flat connector with different color, convenient for clinical use.

Four claws flat connector can hold bone surface, ensure the wire position stable when tightening.

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
■ 9392 扁接头压紧钳
Flat connector crimper

档位调节
Adjust the gears

- 不同规格钛缆选择不同档位（如图）
- 一次性压紧扁接头，无松脱现象
- 钛合金材质，轻便
- 有保险装置，确保扁接头压到位后，保险装置弹开
- Different cable size with different gear (as shown).
- One-time compress the flat connector, no loose phenomenon.
- Made of titanium alloy, handiness.
- With a safety device, ensure to sprang open after apposite compression.

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

Easy to use, one-way tighten cable.

The front and back pressure cable lock can be substituted. Tighten cable tool is universal to all cable sizes.

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Prosthesis and
revision femur locking plate

产品编号 Item No.	规格 Specification		
3757		2孔	125mm
3758	左 L	11孔	270mm
3759	右 R	11孔	270mm

配用 $\Phi 5.0$ 锁定钉、 $\Phi 4.5$ 皮质骨螺钉
Used for $\Phi 5.0$ locking screw, $\Phi 4.5$ cortical screw

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Take prosthesis and revision femur locking plate for example to introduce titanium binding system

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Prosthesis and revision femur locking plate has two specifications --- 2 holes and 11 holes. Generally, 2 holes plate with 4 cables and all cables must pass through the plate's cable holes; when using 11 holes plate, cables on the trochanter of femur need pass through the plate's cable holes, cables on the femur diaphysis need match the R grooves.

Before fixing plate, please connect cable and flat connector together, then pass the cable through plate's cable hole. Otherwise, it maybe hard to through wires owing to soft tissues.

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Lead cable guide against bone surface from back to front(avoid injury of soft tissue and blood vessels). Insert cable into the cable guide head, then withdrawal the guide (picture 3). Forward the cable through another hole of flat connector(flat connector has four claws to hold bone surface). The left cable will be in front (picture 4).

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Φ 1.8 titanium cable (Φ 1.9 flat connector) can pass through the plate's cable hole from both side. When apply to the common femoral lateral approach, recommend passing from the back to front (picture 1). Use beat tool to knock the hook into greater trochanter (picture 2).

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Pass all the left cable from front pressure cable lock to the back one. Straighten cable and lock the later pressure lock, press the trigger (picture 6). Lock the front pressure lock upon getting the expected tension (picture 7). Open the back pressure lock, remove the tighten tool (picture 8). (Remark: don't remove the front pressure lock!)

Pay attention to the angle of tighten tool and flat connector, as parallel as possible. To avoid the cable stuck and guarantee getting expected tension.

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After finishing all cables, owing to the front pressure cable locks remain locked. Please tighten again according to above steps if needed to get satisfying tension.

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Generally, recommend no tighten again (unless cable loose).

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Adjust the gear to $\phi 1.8$ (picture 10), compress all flat connectors one by one(keep compressing over 10 seconds, picture 11). Then open the front pressure cable lock and remove (picture 12).

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Cut titanium cable along the edge of flat connector(please cut one-time, if cut again, must lay gauze to avoid scraps of titanium cable). Then screwing in $\phi 5.0$ locking screws according to specific circumstances (matching $\phi 5.0$ locking instrument set).


图 13 Picture 13

图 14 Picture 14

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INNOVATING CONSTANTLY
KEEP STRIVING FOR EXCELLENCE

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ADD：Shuangyang Building, Baoji, Yangshe Town, Zhangjiagang City, Jiangsu Province, China

电话 Telephone：0512-58268332 传真 FAX：0512-58268337 邮箱 E-Mail：shuangyang@jsshuangyang.com

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Because of the printing, it may vary with the color of the actual product, please subject to the product color when using.

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Only for the medical professional technical personnel and the clinician reference reading, the final interpretation right of this event is reserved by Jiangsu Shuangyang Medical Instrument Co., Ltd.