




LAMINOPLASTY FIXATION SYSTEM

椎板固定板系统



江苏双羊医疗器械有限公司

Jiangsu Shuangyang
Medical Instrument Co., Ltd.



江苏双羊医疗器械有限公司
Jiangsu Shuangyang Medical Instrument Co., Ltd.



企业简介

Company Introduction

江苏双羊医疗器械有限公司，创建于2001年，占地面积18000平方米，建筑面积15000多平方米，注册资金2000万。是一家从事骨科植入物产品的研发、制造、销售和服务的民族企业，获得多项国家级专利。

公司选择钛及钛合金为原材料，严把质量关，选择国内外知名企业宝钛、ZAPP等作为原材料供应商，并配备世界一流的加工中心、纵切车床、数控铣床、超声波清洗等生产设备和装置，以及万能试验机、电子扭转机、数字投影仪等精密的量具。公司拥有完善的管理体系，先后通过了ISO9001:2015质量管理体系认证，ISO13485:2016医疗器械质量管理体系认证，TUV的CE认证，并于2007年率先通过国家局组织的《医疗器械生产质量管理规范植入性医疗器械实施细则（试点）》检查。



公司针对人体各个骨骼部位的特点，在著名骨科专家、教授、临床医师的悉心指导与帮助下，推出众多主导产品：锁定接骨板固定系统、钛质接骨板固定系统、钛质空心接骨螺钉及垫片、钛质胸肋骨系统、锁定颌面内固定系统、颌面内固定系统、钛质捆绑系统、解剖型钛网系统、基础工具系列等，并拥有专业的配套手术器械包，以满足临床不同的需求。公司的产品因设计可靠、加工精细、医生临床使用方便、患者愈合期短而受到了许多临床医师和患者的高度评价。

中国梦，双羊梦！双羊作为一个有使命感、有责任心、有抱负、有人文情怀的公司的初衷不会变，始终以“以人为本，诚信至上，不断创新，追求卓越”的理念不会变，争创医疗器械领域民族品牌领头羊的决心不会变。愿更多有志者与双羊共同创造美好的未来！



Jiangsu Shuangyang Medical Instrument Co., Ltd. was founded in 2001, covers an area of 18,000 m², including a floor area of over 15,000 m². Its registered capital reaches 20 million Yuan. As a national enterprise dedicated to the R&D, manufacturing, sales and service of orthopaedic implants, we have obtained several national patents.


Titanium and titanium alloys are our raw materials. We perform strict quality control, and select domestic and international famous brands, such as Baoti and ZAPP, as our raw material suppliers. Meanwhile, we are equipped with world class production equipment and devices including machining center, slitting lathe, CNC milling machine, and ultrasonic cleaner, etc., as well as precise measuring implements including universal tester, electronic torsion tester and digital projector, etc. Thanks to a sophisticated management system, we have acquired ISO9001: 2015 Certificate of Quality Management System, ISO13485:2016 Certificate of Quality Management System for Medical Devices, and CE certificate of TUV. We are also the first to pass the inspection according to the *Enforcement Regulation (Pilot) for Implantable Medical Devices of Good Manufacturing Practice for Medical Devices* organized by the National Bureau in 2007.

Thanks to meticulous guidance and supports from distinguished orthopaedic specialists, professors and clinicians, we have launched numerous leading products customized for different human skeletal parts, including locking bone plate fixation system, titanium bone plate fixation system, titanium cannulated bone screw & gasket, titanium sternocostal system, locking maxillofacial internal fixation system, maxillofacial internal fixation system, titanium binding system, anatomic titanium mesh system, and basic tool series, etc. We also have professional supporting surgical instrument sets to meet various clinical needs. Extensive acclaims have been received from clinicians and patients for our easy-to-use products with reliable design and fine machining, which can bring a short healing period.

China dream and Shuangyang dream! We will stick to our original intention to be a mission-driven, responsible, ambitious and humanistic company, and adhere to our idea of “people orientation, integrity, innovation, and excellence”. We are determined to be a leading national brand in the medical Instrument industry. At Shuangyang, we always welcome aspiring talents to co-create a bright future with us.

产品特点 Product Features

- 接骨板预裁、预弯、多样化的设计，提供临床更多的选择；
 - 接骨板采用表面阳极氧化工艺，提高产品表面的硬度和耐磨性；
 - 接骨板采用下沉孔设计，板钉结合更紧密，达到低切迹的修复效果，减少软组织不适感；
 - 接骨板边缘设计圆钝，减少对周围软组织的刺激。
-
- The pre-cut, pre-bending and diversified design of plate provide more clinical options.
 - The plate surface adopt anodizing technology, can enhance surface hardness and abrasive resistance.
 - The plate adopt sinking hole design, to make the plate and screw combine closer to achieve low-profile fixation and reduce the discomfort of soft tissue.
 - The edges of plate are rounded to reduce irritation to the surrounding soft tissue.



01

自钻 / 自攻螺钉，十字驱动，持钉更方便
Self drilling/self tapping screw with cross head, more convenient to hold screws.

02


螺钉根部特有的加强设计有效防止断钉。
Root of the screw has special reinforce design to effectively prevent screw broken.

03

抗拔出螺牙角度设计，确保螺钉的把持力。
Screw teeth have anti-pull angle design, ensure the pullout strength of the screw.

04

特有的同心度控制，确保进钉容易，上钉轻松。
Special concentricity control can ensure screwing in easily.



椎板成形板 Laminoplasty system mini plate



直型板 Straight Plate Type I

订货编码 Item No.	规格 Specification	
13.29.01.32006220	32mm	6 孔
13.29.01.64012220	64mm	12 孔



直型板 Straight Plate Type II

订货编码 Item No.	规格 Specification	
13.29.01.55012221	55mm	12 孔



单弯型板 Single bend plate

订货编码 Item No.	规格 Specification	
13.29.03.25004220	25mm	4 孔
13.29.03.27004220	27mm	4 孔
13.29.03.29004220	29mm	4 孔
13.29.03.31004220	31mm	4 孔
13.29.03.33004220	33mm	4 孔



双弯型板 Double Bend Plate Type



订货编码 Item No.	规格 Specification	
13.29.06.06004220	6mm	4 孔
13.29.06.07504220	7.5mm	4 孔
13.29.06.09004220	9mm	4 孔
13.29.06.10504220	10.5mm	4 孔
13.29.06.12004220	12mm	4 孔

纵脚开口型板 Wide Mouth Plate Type



订货编码 Item No.	规格 Specification		
13.29.07.08004221	8mm	4 孔	
13.29.07.10004221	10mm	4 孔	
13.29.07.11504221	11.5mm	4 孔	
13.29.07.08004222	8mm	4 孔	6.8
13.29.07.10004222	10mm	4 孔	
13.29.07.11504222	11.5mm	4 孔	

横脚开口型板 Open Door Plate Type



订货编码 Item No.	规格 Specification		
13.29.08.07504221	7.5mm	4 孔	
13.29.08.09504221	9.5mm	4 孔	
13.29.08.11504221	11.5mm	4 孔	
13.29.08.07504222	7.5mm	4 孔	6.8
13.29.08.09504222	9.5mm	4 孔	
13.29.08.11504222	11.5mm	4 孔	


椎板螺钉 Laminoplasty system screw

自攻螺钉
2.0 Self-tapping
screw




订货编码 Item No.	规格 Specification
13.29.10.20050110	Ø 2.0×5mm
13.29.10.20070110	Ø 2.0×7mm
13.29.10.20090110	Ø 2.0×9mm
13.29.10.20110110	Ø 2.0×11mm

自攻螺钉
2.4 Self-tapping
screw



订货编码 Item No.	规格 Specification
13.29.10.24060110	Ø 2.4×6mm
13.29.10.24080110	Ø 2.4×8mm
13.29.10.24100110	Ø 2.4×10mm
13.29.10.24120110	Ø 2.4×12mm

自攻自钻螺钉
2.0 Self-drilling
screw



订货编码 Item No.	规格 Specification
13.29.20.20050110	Ø 2.0×5mm
13.29.20.20070110	Ø 2.0×7mm
13.29.20.20090110	Ø 2.0×9mm
13.29.20.20110110	Ø 2.0×11mm



不断创新 · 追求卓越
INNOVATING CONSTANTLY
KEEP STRIVING FOR EXCELLENCE

www.jsshuangyang.com

地址：江苏省张家港市杨舍镇包基双羊大厦 邮编 Postcode: 215600

ADD: Shuangyang Building, Baoji, Yangshe Town, Zhangjiagang City, Jiangsu Province, China

电话 Telephone: 0512-58268332 传真 FAX: 0512-58268337 E-Mail: shuangyang@jsshuangyang.com



微信号：shuangyang8

由于印刷的关系，可能与实际产品的颜色有所差异，使用时请以实物颜色为准。

Because of the printing, it may vary with the color of the actual product, please subject to the product color when using.

仅供医疗器械专业技术人员、临床医师参考阅读，最终解释权归江苏双羊医疗器械有限公司所有。

Only for the medical professional technical personnel and the clinician reference reading, the final interpretation right of this event is reserved by Jiangsu Shuangyang Medical Instrument Co., Ltd.